**IDROSSILAPATITE MICROMETRICA AT®****A RAPIDO INDURIMENTO**

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**IDROSSILAPATITE MICROMETRICA AT®****FAST HARDENING**

Granular osteoconductive micrometric hydroxylapatite with fast hardening osteoconductive for guided tissue regeneration in dental surgery

"Guided regeneration materials system  
Collagene AT® Idrossilapatite AT®"

**Properties and Technical Features**

The use of fast hardening **Idrossilapatite micrometrica AT®** increases the primary retention of implants and prevents epithelial and connective tissue cells from interfering with bone growth and periodontal tissue growth. Surface micro-irregularity of **Idrossilapatite micrometrica AT®** allows the absorption of plasmatic proteins with function of growth factors and with osteoinductive action.

**Idrossilapatite micrometrica AT®**, calcium hydroxylapatite with molecular weight 502,32 has been designed and created because its particular granulometry with the addition of calcium sulphate and calcium phosphate allows a better manipulation and application of the material into the operating field, thanks to its reciprocal adhesion of the particles.

Sphere shaped granulometry adhere each other for capillarity and leave continuous intermediate spaces, which have an osteoconductive action and are progressively colonized by bony trabeculations.

Reabsorbable progressively in centripetal direction, towards the inside of the granule in a time span of 4-10 weeks, during the formation of new bone.

**Main Utilization**

- ♦ To increase the primary retention of implants.

**Modalities**

- 1) pour **Idrossilapatite micrometrica AT®** in a sterile glass and add physiologic solution stirring until the surface of the mixture appears just glossy
- 2) apply **Idrossilapatite micrometrica AT®** on the lateral surface and on the extremity of the implant
- 3) screw immediately the implant and remove the hydroxylapatite in excess.

- ♦ The **Idrossilapatite micrometrica AT®** rapid curing has a 1-3 minutes of processing time and a consolidation time of 4-45 minutes

- ♦ The fast hardening **Idrossilapatite micrometrica AT®** on the internal part of thread allows to obtain the primary stability of implants during the first osteointegration, which occurs in correspondence to the external part of thread that remains in direct contact with the bone tissue. When the micrometric hydroxylapatite is reabsorbed after 4-10 weeks the osteointegration extends in the internal zone of the implants thread.

- ♦ The fast hardening **Idrossilapatite micrometrica AT®** is used to increase the torque of insertion of implants up to 100%

- ♦ The fast hardening **Idrossilapatite micrometrica AT®** is the solution to obtain the immediate optimal stability of the implants in the bone tissue type 3 and 4.

**Packaging and Conservation**

The fast hardening **Idrossilapatite micrometrica AT®** is sold in:

- ♦ 2 glass bottles of 0,5g single use (**code 2607**) with bakelite cap, and blister with double welding.

The fast hardening **Idrossilapatite micrometrica AT®** is sterilized by gamma radiation.

**Precautions**

The **Idrossilapatite micrometrica AT®** is contraindicated in patients with acute infections or contaminated wounds.

Medical device **Idrossilapatite micrometrica AT®** must be used only by Dental Surgeons, Physicians, Dentists qualified in surgical dental procedures.

This medical device can be used only if the packaging is intact.

Adhesive labels included for traceability to be affixed in the patient's file.

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