# COLLAGENE AT®

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# COLLAGENE AT®

Resorbable Collagen Membrane for guided tissue regeneration in dental surgery

“Guided regeneration materials system  
Collagene AT® Idrossilapatite AT®”

## Properties and technical features

Resorbable membrane **Collagene AT®** is indicated in every case where guided tissue regeneration is needed, in Periodontology, in Implantology, in Oral Surgery and for the prevention of post-surgical reabsorption of the alveolar ridges.

**Collagene AT®** membrane is made of pure, and lyophilized collagen of equine extraction.

**Collagene AT®** is a membrane with micro-irregularities at the surface that creates a barrier effect; amongst the collagen layers with reticular links there are micro spaces of 10-20 microns avoiding cells migration.

The barrier effect prevents epithelial and soft connective tissue cells from interfering with bone growth or paradontal tissue growth. The maintenance of a space beneath the membrane allows osteoblasts migration and new bone growth or new paradontal growth.

Membranes create a barrier effect for the first six weeks.

Membranes are totally reabsorbed within 180 days and therefore eliminating the need for second surgery.

For an optimal guided tissue regeneration, in order to maintain space within the defect cavity for bone regeneration, a grafting material should always be used like **Idrossilapatite AT®**.

**Collagene AT®** membrane is compatible with all the other different bone substitutes available in the market.

## Packaging and conservation

Each box contains 6 **Collagene AT®** membranes, packaged individually in transparent 0,8 mm thick PET G blister, flexible and resistant, unshrinkable at high temperatures and sealed with a TYVEC lid in order to ensure the absolute sterility of the product.

They are sterilized by gamma radiation.

**Collagene AT®** membranes are for single use and must be stored at temperature below 40°C. Do not sterilize in autoclave.

## Handling

Sterile instruments are needed.

Open the blister containing the membrane removing the Tyvec lid over a sterile sheet.

To obtain a better result, the membrane has to overlap the bony defect for 2-3 mm in healthy bone.

At the beginning the membrane looks rigid but as soon as hydrated with physiological solution or blood, it becomes manageable without losing its form and texture.

Adapt the membrane with dry sterile gauzes.

In case of membrane exposure, it is recommended to protect it with a gel containing clorexidine and surgical dressing until the epithelium covers the area.

**Collagene AT®** can be used in filling of postextraction sockets, in filling of paradontal pockets, in maxillary sinus lifts, in sinus-oral communications, in atrophic ridge reconstruction, in periapical surgery after apicectomies and in every case of oral surgery where bone reconstruction is needed.

## Precautions

**Collagene AT®** is contraindicated in patients with acute infections, contaminated wounds or intolerance to equine tissues.

This medical device can be used only by Dental Surgeons, Physicians, Dentists qualified in surgical dental procedures.

Before using the **Collagene AT®**, make sure the thermo sealed blister is intact.

Adhesive labels included for traceability to be affixed in the patient's file.

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